

NEXTflex[™] ChIP-Seq Barcode Indices

NEXTflex [™] ChIP-Seq Barcodes	Index Sequence
ChIP-Seq Adapter 1	CGATGT
ChIP-Seq Adapter 2	TGACCA
ChIP-Seq Adapter 3	ACAGTG
ChIP-Seq Adapter 4	GCCAAT
ChIP-Seq Adapter 5	CAGATC
ChIP-Seq Adapter 6	CTTGTA
ChIP-Seq Adapter 7	ATCACG
ChIP-Seq Adapter 8	TTAGGC
ChIP-Seq Adapter 9	ACTTGA
ChIP-Seq Adapter 10	GATCAG
ChIP-Seq Adapter 11	TAGCTT
ChIP-Seq Adapter 12	GGCTAC
ChIP-Seq Adapter 13	AGTCAA
ChIP-Seq Adapter 14	AGTTCC
ChIP-Seq Adapter 15	ATGTCA
ChIP-Seq Adapter 16	CCGTCC
ChIP-Seq Adapter 17	GTAGAG
ChIP-Seq Adapter 18	GTCCGC
ChIP-Seq Adapter 19	GTGAAA
ChIP-Seq Adapter 20	GTGGCC
ChIP-Seq Adapter 21	GTTTCG
ChIP-Seq Adapter 22	CGTACG
ChIP-Seq Adapter 23	GAGTGG
ChIP-Seq Adapter 24	GGTAGC
ChIP-Seq Adapter 25	ACTGAT
ChIP-Seq Adapter 26	ATGAGC
ChIP-Seq Adapter 27	ATTCCT
ChIP-Seq Adapter 28	CAAAAG
ChIP-Seq Adapter 29	CAACTA
ChIP-Seq Adapter 30	CACCGG
ChIP-Seq Adapter 31	CACGAT
ChIP-Seq Adapter 32	CACTCA
ChIP-Seq Adapter 33	CAGGCG
ChIP-Seq Adapter 34	CATGGC
ChIP-Seq Adapter 35	CATTTT
ChIP-Seq Adapter 36	CCAACA
ChIP-Seq Adapter 37	CGGAAT
ChIP-Seq Adapter 38	CTAGCT
ChIP-Seq Adapter 39	CTATAC
ChIP-Seq Adapter 40	CTCAGA

ChIP-Seq Adapter 41	GCGCTA
ChIP-Seq Adapter 42	TAATCG
ChIP-Seq Adapter 43	TACAGC
ChIP-Seq Adapter 44	TATAAT
ChIP-Seq Adapter 45	TCATTC
ChIP-Seq Adapter 46	TCCCGA
ChIP-Seq Adapter 47	TCGAAG
ChIP-Seq Adapter 48	TCGGCA

©Bio Scientific Corp. • 2015 v1.1