

## NEXTflex-96™ ChIP-Seq Barcode Indices

NEXTflex-96™ ChIP-Seq Barcodes	Index Sequence
NEXTflex-96™ ChIP-Seq Adapter 1	AACGTGAT
NEXTflex-96™ ChIP-Seq Adapter 2	AAACATCG
NEXTflex-96™ ChIP-Seq Adapter 3	ATGCCTAA
NEXTflex-96™ ChIP-Seq Adapter 4	AGTGGTCA
NEXTflex-96™ ChIP-Seq Adapter 5	ACCACTGT
NEXTflex-96™ ChIP-Seq Adapter 6	ACATTGGC
NEXTflex-96™ ChIP-Seq Adapter 7	CAGATCTG
NEXTflex-96™ ChIP-Seq Adapter 8	CATCAAGT
NEXTflex-96™ ChIP-Seq Adapter 9	CGCTGATC
NEXTflex-96™ ChIP-Seq Adapter 10	ACAAGCTA
NEXTflex-96™ ChIP-Seq Adapter 11	CTGTAGCC
NEXTflex-96™ ChIP-Seq Adapter 12	AGTACAAG
NEXTflex-96™ ChIP-Seq Adapter 13	ACAACCA
NEXTflex-96™ ChIP-Seq Adapter 14	AACCGAGA
NEXTflex-96™ ChIP-Seq Adapter 15	AACGCTTA
NEXTflex-96™ ChIP-Seq Adapter 16	AAGACGGA
NEXTflex-96™ ChIP-Seq Adapter 17	AAGGTACA
NEXTflex-96™ ChIP-Seq Adapter 18	ACACAGAA
NEXTflex-96™ ChIP-Seq Adapter 19	ACAGCAGA
NEXTflex-96™ ChIP-Seq Adapter 20	ACCTCAA
NEXTflex-96™ ChIP-Seq Adapter 21	ACGCTCGA
NEXTflex-96™ ChIP-Seq Adapter 22	ACGTATCA
NEXTflex-96™ ChIP-Seq Adapter 23	ACTATGCA
NEXTflex-96™ ChIP-Seq Adapter 24	AGAGTCAA
NEXTflex-96™ ChIP-Seq Adapter 25	AGATCGCA
NEXTflex-96™ ChIP-Seq Adapter 26	AGCAGGAA
NEXTflex-96™ ChIP-Seq Adapter 27	AGTCACTA
NEXTflex-96™ ChIP-Seq Adapter 28	ATCCTGTA
NEXTflex-96™ ChIP-Seq Adapter 29	ATTGAGGA
NEXTflex-96™ ChIP-Seq Adapter 30	CAACCACA
NEXTflex-96™ ChIP-Seq Adapter 31	GACTAGTA
NEXTflex-96™ ChIP-Seq Adapter 32	CAATGGAA
NEXTflex-96™ ChIP-Seq Adapter 33	CACTTCGA
NEXTflex-96™ ChIP-Seq Adapter 34	CAGCGTTA
NEXTflex-96™ ChIP-Seq Adapter 35	CATACCAA
NEXTflex-96™ ChIP-Seq Adapter 36	CCAGTTCA
NEXTflex-96™ ChIP-Seq Adapter 37	CCGAAGTA
NEXTflex-96™ ChIP-Seq Adapter 38	CCGTGAGA
NEXTflex-96™ ChIP-Seq Adapter 39	CCTCCTGA
NEXTflex-96™ ChIP-Seq Adapter 40	CGAACTTA

NEXTflex-96™ ChIP-Seq Adapter 41	CGACTGGA
NEXTflex-96™ ChIP-Seq Adapter 42	CGCATACA
NEXTflex-96™ ChIP-Seq Adapter 43	CTCAATGA
NEXTflex-96™ ChIP-Seq Adapter 44	CTGAGCCA
NEXTflex-96™ ChIP-Seq Adapter 45	CTGGCATA
NEXTflex-96™ ChIP-Seq Adapter 46	GAATCTGA
NEXTflex-96™ ChIP-Seq Adapter 47	CAAGACTA
NEXTflex-96™ ChIP-Seq Adapter 48	GAGCTGAA
NEXTflex-96™ ChIP-Seq Adapter 49	GATAGACA
NEXTflex-96™ ChIP-Seq Adapter 50	GCCACATA
NEXTflex-96™ ChIP-Seq Adapter 51	GCGAGTAA
NEXTflex-96™ ChIP-Seq Adapter 52	GCTAACGA
NEXTflex-96™ ChIP-Seq Adapter 53	GCTCGGTA
NEXTflex-96™ ChIP-Seq Adapter 54	GGAGAACA
NEXTflex-96™ ChIP-Seq Adapter 55	GGTGCGAA
NEXTflex-96™ ChIP-Seq Adapter 56	GTACGCAA
NEXTflex-96™ ChIP-Seq Adapter 57	GTCGTAGA
NEXTflex-96™ ChIP-Seq Adapter 58	GTCTGTCA
NEXTflex-96™ ChIP-Seq Adapter 59	GTGTTCTA
NEXTflex-96™ ChIP-Seq Adapter 60	TAGGATGA
NEXTflex-96™ ChIP-Seq Adapter 61	TATCAGCA
NEXTflex-96™ ChIP-Seq Adapter 62	TCCGTCTA
NEXTflex-96™ ChIP-Seq Adapter 63	TCTTCACA
NEXTflex-96™ ChIP-Seq Adapter 64	TGAAGAGA
NEXTflex-96™ ChIP-Seq Adapter 65	TGGAACAA
NEXTflex-96™ ChIP-Seq Adapter 66	TGGCTTCA
NEXTflex-96™ ChIP-Seq Adapter 67	TGGTGGTA
NEXTflex-96™ ChIP-Seq Adapter 68	TTCACGCA
NEXTflex-96™ ChIP-Seq Adapter 69	AACTCACC
NEXTflex-96™ ChIP-Seq Adapter 70	AAGAGATC
NEXTflex-96™ ChIP-Seq Adapter 71	AAGGACAC
NEXTflex-96™ ChIP-Seq Adapter 72	AATCCGTC
NEXTflex-96™ ChIP-Seq Adapter 73	AATGTTGC
NEXTflex-96™ ChIP-Seq Adapter 74	ACACGACC
NEXTflex-96™ ChIP-Seq Adapter 75	ACAGATTC
NEXTflex-96™ ChIP-Seq Adapter 76	AGATGTAC
NEXTflex-96™ ChIP-Seq Adapter 77	AGCACCTC
NEXTflex-96™ ChIP-Seq Adapter 78	AGCCATGC
NEXTflex-96™ ChIP-Seq Adapter 79	AGGCTAAC
NEXTflex-96™ ChIP-Seq Adapter 80	ATAGCGAC
NEXTflex-96™ ChIP-Seq Adapter 81	ATCATTCC
NEXTflex-96™ ChIP-Seq Adapter 82	ATTGGCTC
NEXTflex-96™ ChIP-Seq Adapter 83	CAAGGAGC
NEXTflex-96™ ChIP-Seq Adapter 84	CACCTTAC
NEXTflex-96™ ChIP-Seq Adapter 85	CCATCCTC
NEXTflex-96™ ChIP-Seq Adapter 86	CCGACAAC
NEXTflex-96™ ChIP-Seq Adapter 87	CCTAATCC

NEXTflex-96™ ChIP-Seq Adapter 88	CCTCTATC
NEXTflex-96™ ChIP-Seq Adapter 89	CGACACAC
NEXTflex-96™ ChIP-Seq Adapter 90	CGGATTGC
NEXTflex-96™ ChIP-Seq Adapter 91	CTAAGGTC
NEXTflex-96™ ChIP-Seq Adapter 92	GAACAGGC
NEXTflex-96™ ChIP-Seq Adapter 93	GACAGTGC
NEXTflex-96™ ChIP-Seq Adapter 94	GAGTTAGC
NEXTflex-96™ ChIP-Seq Adapter 95	GATGAATC
NEXTflex-96™ ChIP-Seq Adapter 96	GCCAAGAC