

NEXTflex-HT[™] Barcode Indices

NEXTflex-HT [™] Barcodes	Sequence (5' - 3')
Barcode Adapter 1	ATCTAGCCGGCC
Barcode Adapter 2	TATCTCTTCCTT
Barcode Adapter 3	TAGATGCCGTCC
Barcode Adapter 4	CGCTCCTTCCTT
Barcode Adapter 5	CGAGCGCCGTCC
Barcode Adapter 6	TCGAGAGTCCGG
Barcode Adapter 7	GCTCGGCCGGCC
Barcode Adapter 8	AGCTCAGTCCGG
Barcode Adapter 9	ACTGCCGCGTCC
Barcode Adapter 10	TTACTTCTCCGG
Barcode Adapter 11	TGACGCGCGGCC
Barcode Adapter 12	AATGATCTCCGG
Barcode Adapter 13	AAGCGTACGTCC
Barcode Adapter 14	TTATCAGTCCTT
Barcode Adapter 15	GTCATCGCGTCC
Barcode Adapter 16	CCGTCTCTCCGG
Barcode Adapter 17	TTCGCTACGTCC
Barcode Adapter 18	ACGCTCTTCCGG
Barcode Adapter 19	CAGTACGCGGCC
Barcode Adapter 20	GTACTGATCCTT
Barcode Adapter 21	GGATATACGGCC
Barcode Adapter 22	CATAGCTTCCGG
Barcode Adapter 23	GAGGCATCGGCC
Barcode Adapter 24	ATAAGGATCCGG
Barcode Adapter 25	CTCCGATCGTCC
Barcode Adapter 26	TATTCGATCCGG
Barcode Adapter 27	TCTTAATCGTCC
Barcode Adapter 28	CGCCTGATCCGG
Barcode Adapter 29	AGAATATCGTCC
Barcode Adapter 30	GCGGAGATCCGG
Barcode Adapter 31	CATGGCCGGTCC
Barcode Adapter 32	GTCACGTACCTT
Barcode Adapter 33	TGCAACCGGTCC
Barcode Adapter 34	CCTGTTGACCGG
Barcode Adapter 35	AGCCCAAGGTCC
Barcode Adapter 36	GCTTGGTACCGG
Barcode Adapter 37	TCGGGAAGGTCC
Barcode Adapter 38	CGAACGTACCGG
Barcode Adapter 39	GATTTAAGGGCC
Barcode Adapter 40	AGAACTGACCTT

Barcode Adapter 41	CCGCCTTGGTCC
Barcode Adapter 42	TGATGCAACCGG
Barcode Adapter 43	GGCGGTTGGGCC
Barcode Adapter 44	ACTACCAACCGG
Barcode Adapter 45	AATAATTGGTCC
Barcode Adapter 46	GTCGTCAACCGG
Barcode Adapter 47	AAATCTCAGGCC
Barcode Adapter 48	GTGCGCGGCCGG
Barcode Adapter 49	TTAGTCAGTCC
Barcode Adapter 50	CACGCCGGCCGG
Barcode Adapter 51	CCCGATCAGTCC
Barcode Adapter 52	TGTATCGGCCGG
Barcode Adapter 53	CTTTCAGAGTCC
Barcode Adapter 54	TACCGGCGGCCGG
Barcode Adapter 55	GAAAGAGAGTCC
Barcode Adapter 56	CTGGCTAGCCTT
Barcode Adapter 57	AGGGAAGAGGCC
Barcode Adapter 58	GACCGTAGCCTT
Barcode Adapter 59	TCCCTAGAGTCC
Barcode Adapter 60	AGTTATAGCCTT
Barcode Adapter 61	GCCTCGAAGTCC
Barcode Adapter 62	CGTCGCGGCCTT
Barcode Adapter 63	CGGAGGAAGGCC
Barcode Adapter 64	GCAGCCGGCCTT
Barcode Adapter 65	TAAGAGAAGTCC
Barcode Adapter 66	ATGATCGGCCTT
Barcode Adapter 67	ATTCTGAAGTCC
Barcode Adapter 68	TACTACGGCCTT
Barcode Adapter 69	TGGTCCTAGTCC
Barcode Adapter 70	AACATTAGCCGG
Barcode Adapter 71	ACCAGCTAGTCC
Barcode Adapter 72	TTGTATAGCCGG
Barcode Adapter 73	GTTGACTAGTCC
Barcode Adapter 74	CCACGTAGCCGG
Barcode Adapter 75	CAACTCTAGTCC
Barcode Adapter 76	TGTGCGCGCCTT
Barcode Adapter 77	GGTACTGTGGCC
Barcode Adapter 78	AAATTAACCCTT
Barcode Adapter 79	TTGCATGTGGCC
Barcode Adapter 80	CCCGGAACCCTT
Barcode Adapter 81	AACGTTGTGTCC
Barcode Adapter 82	GGGCCAACCCCTT
Barcode Adapter 83	CACACCATGTCC
Barcode Adapter 84	TTGGTTCCCGG
Barcode Adapter 85	GTGTGCATGTCC
Barcode Adapter 86	CCCAATTCCCGG
Barcode Adapter 87	TGTGTCATGTCC

Barcode Adapter 88	AAACCTCCCCGG
Barcode Adapter 89	ATGACGTTGTCC
Barcode Adapter 90	TCCTTAACCCGG
Barcode Adapter 91	CGTCAGTTGTCC
Barcode Adapter 92	GAAGGAACCCGG
Barcode Adapter 93	GCAGTGTTGGCC
Barcode Adapter 94	CTCCAACCCGG
Barcode Adapter 95	GTATGACCAGCC
Barcode Adapter 96	CCTAAGGTTCCGG
Barcode Adapter 97	TGCGTACCATCC
Barcode Adapter 98	AAGCCGGTTCGG
Barcode Adapter 99	ATAACTGCATCC
Barcode Adapter 100	GAGGGCCTTCGG
Barcode Adapter 101	TATTGTGCAGCC
Barcode Adapter 102	AGAAACCTTCGG
Barcode Adapter 103	CGCCATGCATCC
Barcode Adapter 104	TCTTTCCTTCGG
Barcode Adapter 105	GCGGTTGCATCC
Barcode Adapter 106	ATCCCAATTCTT
Barcode Adapter 107	TCGACCACAGCC
Barcode Adapter 108	CTCTTGGTTCTT
Barcode Adapter 109	AGCTGCACAGCC
Barcode Adapter 110	TCTCCGGTTCTT
Barcode Adapter 111	CTAGTCACAGCC
Barcode Adapter 112	GAGAAGGTTCTT
Barcode Adapter 113	CCGTGGTCATCC
Barcode Adapter 114	TGACCAATTCGG
Barcode Adapter 115	TTACAGTCATCC
Barcode Adapter 116	CAGTTAATTCGG
Barcode Adapter 117	AATGTGTCAGCC
Barcode Adapter 118	TTCAACCTTCTT
Barcode Adapter 119	GCTTCTCGATCC
Barcode Adapter 120	ATAATATATCTT
Barcode Adapter 121	CGAAGTCGATCC
Barcode Adapter 122	TCGGCCGATCGG
Barcode Adapter 123	TAGGATCGATCC
Barcode Adapter 124	CTAATCGATCGG
Barcode Adapter 125	ATCCTTCGATCC
Barcode Adapter 126	GATTACGATCGG
Barcode Adapter 127	TGATCAGGAGCC
Barcode Adapter 128	CCGCGGCATCGG
Barcode Adapter 129	ACTAGAGGATCC
Barcode Adapter 130	GTATATAATCTT
Barcode Adapter 131	GTCGAAGGAGCC
Barcode Adapter 132	AATATGCATCGG
Barcode Adapter 133	CAGCTAGGATCC
Barcode Adapter 134	TGCGCTAATCTT

Barcode Adapter 135	TTCAGGAGAGCC
Barcode Adapter 136	ACGTAATATCGG
Barcode Adapter 137	CCTGAGAGATCC
Barcode Adapter 138	TTACGCGATCTT
Barcode Adapter 139	GGA CTGAGATCC
Barcode Adapter 140	CATGCATATCGG
Barcode Adapter 141	CTCTCCTGATCC
Barcode Adapter 142	GATCGGCATCTT
Barcode Adapter 143	GAGAGCTGAGCC
Barcode Adapter 144	ATAGCTAATCGG
Barcode Adapter 145	AGAGACTGATCC
Barcode Adapter 146	TCGATGCATCTT
Barcode Adapter 147	TCTCTCTGAGCC
Barcode Adapter 148	CGCTATAATCGG
Barcode Adapter 149	TCAGGCCAAGCC
Barcode Adapter 150	CGGACTGGTCGG
Barcode Adapter 151	CTGAACCAATCC
Barcode Adapter 152	GCCTGTGGTCGG
Barcode Adapter 153	GACTTCCAAGCC
Barcode Adapter 154	CTTCAGTGTCTT
Barcode Adapter 155	CCACCGGAATCC
Barcode Adapter 156	TGGTGACGTCGG
Barcode Adapter 157	AACAAGGAAGCC
Barcode Adapter 158	GTTGTACGTCGG
Barcode Adapter 159	TTGTTGGAATCC
Barcode Adapter 160	AAACACAGTCTT
Barcode Adapter 161	ATGGGTTAATCC
Barcode Adapter 162	GAAACCAGTCGG
Barcode Adapter 163	GCAAATTAACCC
Barcode Adapter 164	CTTTGCAGTGGG
Barcode Adapter 165	TTTGCGCTATCC
Barcode Adapter 166	CCACTCTCTCTT
Barcode Adapter 167	AAACGGCTATCC
Barcode Adapter 168	TGTGAAGCTCGG
Barcode Adapter 169	GGGTAGCTATCC
Barcode Adapter 170	AACAGCTCTCTT
Barcode Adapter 171	CCCATGCTAGCC
Barcode Adapter 172	TTGTCCTCTCTT
Barcode Adapter 173	GAAGCCGTATCC
Barcode Adapter 174	CTGAGGACTCTT
Barcode Adapter 175	CTTCGCGTAGCC
Barcode Adapter 176	GACTCGACTCTT
Barcode Adapter 177	TCCTACGTATCC
Barcode Adapter 178	AGTCTGACTCTT
Barcode Adapter 179	ATGGCGCGCTGG
Barcode Adapter 180	GCCCTCTAGCAA
Barcode Adapter 181	CGTTAGCGCGGG

Barcode Adapter 182	TAAAGCTAGCAA
Barcode Adapter 183	GCAATGCGCTGG
Barcode Adapter 184	ATTCCTAGCAA
Barcode Adapter 185	GGTGCTAGCTGG
Barcode Adapter 186	CCCAGAGAGCAA
Barcode Adapter 187	TCAGCATGCTGG
Barcode Adapter 188	AGGAGTCAGCAA
Barcode Adapter 189	AGTCGATGCTGG
Barcode Adapter 190	GAAGATCAGCAA
Barcode Adapter 191	GACTAATGCTGG
Barcode Adapter 192	CTTCTTCAGCAA
Barcode Adapter 193	CTGATATGCTGG
Barcode Adapter 194	TCCTCTCAGCAA
Barcode Adapter 195	GAGACACACGGG
Barcode Adapter 196	CTAGGTTGGCAA
Barcode Adapter 197	CTCTGACACTGG
Barcode Adapter 198	TCGAATTGGCAA
Barcode Adapter 199	AGAGTACACTGG
Barcode Adapter 200	GATCCTTGGCAA
Barcode Adapter 201	TTCACTGACGGG
Barcode Adapter 202	CCGTTAAGGCAA
Barcode Adapter 203	AAGTGTGACTGG
Barcode Adapter 204	TTACCAAGGCAA
Barcode Adapter 205	ATCCAGTACGGG
Barcode Adapter 206	TATTTCCGGCAA
Barcode Adapter 207	TAGGTGTACCGG
Barcode Adapter 208	ATAAACCGGTAA
Barcode Adapter 209	CTAGAAGTCCGG
Barcode Adapter 210	GAGATTACGTAA
Barcode Adapter 211	GATCTAGTCCGG
Barcode Adapter 212	AGAGCTACGGAA
Barcode Adapter 213	TTGACCATGTGG
Barcode Adapter 214	CCCTTGCCCAA
Barcode Adapter 215	GCATGACTGTGG
Barcode Adapter 216	ATTAATTCCCAA
Barcode Adapter 217	CGTACACTGTGG
Barcode Adapter 218	GCCGTTCCCAA
Barcode Adapter 219	TTTCTCTAGTGG
Barcode Adapter 220	CCAGCGCGCCAA
Barcode Adapter 221	AAAGACTAGTGG
Barcode Adapter 222	TTGATGCGCCAA
Barcode Adapter 223	GAACTGAAGTGG
Barcode Adapter 224	AGTGCCGGCCAA
Barcode Adapter 225	CTTGAGAAGTGG
Barcode Adapter 226	TCACGCGGCCAA
Barcode Adapter 227	TCCAGGAAGTGG
Barcode Adapter 228	CTGTACGGCCAA

Barcode Adapter 229	AGGTCGAAGTGG
Barcode Adapter 230	GACATCGGCCAA
Barcode Adapter 231	CGGCTAGAGTGG
Barcode Adapter 232	TACGCTAGCCAA
Barcode Adapter 233	GCCGAAGAGTGG
Barcode Adapter 234	CGTATTAGCCAA
Barcode Adapter 235	ATTAGAGAGTGG
Barcode Adapter 236	GCATATAGCCAA
Barcode Adapter 237	TAATCAGAGTGG
Barcode Adapter 238	ATGCGTAGCCAA
Barcode Adapter 239	ACCCTTCAGTGG
Barcode Adapter 240	GTGGCATGCCAA
Barcode Adapter 241	TGGGATCAGTGG
Barcode Adapter 242	CACCGATGCCAA
Barcode Adapter 243	CAAAGTCAGTGG
Barcode Adapter 244	TGTTAATGCCAA
Barcode Adapter 245	GTTTCTCAGTGG
Barcode Adapter 246	ACAATATGCCAA
Barcode Adapter 247	GTAAATTGGTGG
Barcode Adapter 248	ACTTGACACCAA
Barcode Adapter 249	ACGGGTTGGTGG
Barcode Adapter 250	TGAACACACCAA
Barcode Adapter 251	TGCCCTTGGTGG
Barcode Adapter 252	CAGGTACACCAA
Barcode Adapter 253	ATATTAAGGTGG
Barcode Adapter 254	GCTACTGACCAA
Barcode Adapter 255	CGCGGAAGGTGG
Barcode Adapter 256	TAGCATGACCAA
Barcode Adapter 257	GCGCCAAGGTGG
Barcode Adapter 258	CGATGTGACCAA
Barcode Adapter 259	GGCTTCCGGTGG
Barcode Adapter 260	CCTCAGTACCAA
Barcode Adapter 261	TTAGGCCGGTGG
Barcode Adapter 262	AAGACGTACCAA
Barcode Adapter 263	GCTATATCGTGG
Barcode Adapter 264	ATATCTCTCCAA
Barcode Adapter 265	CGATAATCGTGG
Barcode Adapter 266	TATAGTCTCCAA
Barcode Adapter 267	TAGCGATCGTGG
Barcode Adapter 268	CGCGATCTCCAA
Barcode Adapter 269	ATCGCATCGTGG
Barcode Adapter 270	GCGCTTCTCCAA
Barcode Adapter 271	TGAATTACGCGG
Barcode Adapter 272	CATTCAGTCTAA
Barcode Adapter 273	GGAAGAGGACGG
Barcode Adapter 274	ACGGCGCATTCC
Barcode Adapter 275	CCTTCAGGAGGG

Barcode Adapter 276	GTAATGCATCCC
Barcode Adapter 277	GAGCTTCGACGG
Barcode Adapter 278	CTATAATATTA
Barcode Adapter 279	CTCGATCGACGG
Barcode Adapter 280	TATATCGATGCC
Barcode Adapter 281	TCTAGTCGATGG
Barcode Adapter 282	CGCGCCGATCCC
Barcode Adapter 283	AGATCTCGACGG
Barcode Adapter 284	GCGCGCGATGCC
Barcode Adapter 285	GTAGTGTCATGG
Barcode Adapter 286	CCTCCAATTCCC
Barcode Adapter 287	CATCAGTCAGGG
Barcode Adapter 288	GTCTTCCTCAA
Barcode Adapter 289	ACGACGTCACGG
Barcode Adapter 290	TTCTTAATTGCC
Barcode Adapter 291	TATGTCACATGG
Barcode Adapter 292	ATCAAGGTTCAA
Barcode Adapter 293	ATACACACATGG
Barcode Adapter 294	TAGTTGGTTCAA
Barcode Adapter 295	AGCGTTGCACGG
Barcode Adapter 296	GAGCCAATTGAA
Barcode Adapter 297	TCGCATGCATGG
Barcode Adapter 298	CTCGGAATTCAA
Barcode Adapter 299	CTATGTGCACGG
Barcode Adapter 300	GCTAACCTTTCC
Barcode Adapter 301	GATACTGCACGG
Barcode Adapter 302	ATCGGCCTTTCC
Barcode Adapter 303	CCGGTACCATGG
Barcode Adapter 304	GTCCCGGTTCCC
Barcode Adapter 305	AATTGACCACGG
Barcode Adapter 306	TGAAAGGTTTCC
Barcode Adapter 307	CCAGTCATGCGG
Barcode Adapter 308	TTCCGGCCGAA
Barcode Adapter 309	CGGCCTTCCCAA
Barcode Adapter 310	GACGTCATGTTT
Barcode Adapter 311	ACACACATGTCC
Barcode Adapter 312	TGGTTGGCCCTT
Barcode Adapter 313	GCTAATCTGGAA
Barcode Adapter 314	CTATGCGCCATT
Barcode Adapter 315	CGATTTCTGAAA
Barcode Adapter 316	GCGCAATCCGGG
Barcode Adapter 317	CAGGGAGTGGAA
Barcode Adapter 318	TGCCATACCAGG
Barcode Adapter 319	AAGAAGATGGAA
Barcode Adapter 320	GTAGTATCCATT
Barcode Adapter 321	TTCTTGATGAAA
Barcode Adapter 322	ACGACATCCGTT

Barcode Adapter 323	AGACCCTTGAAA
Barcode Adapter 324	TATGTTACCGTT
Barcode Adapter 325	TCTGGCTTGAAA
Barcode Adapter 326	CGCACTACCGTT
Barcode Adapter 327	GAGTTCTTGGAA
Barcode Adapter 328	ATACATACCATT
Barcode Adapter 329	AGGCCAGCAGAA
Barcode Adapter 330	TACGTGCTTATT
Barcode Adapter 331	TCCGGAGCAGAA
Barcode Adapter 332	ATGCAGCTTATT
Barcode Adapter 333	GAATTAGCAGAA
Barcode Adapter 334	CGTACGCTTATT
Barcode Adapter 335	TAACCGACATAA
Barcode Adapter 336	AGTGTATTTATT
Barcode Adapter 337	ATTGGGACAGAA
Barcode Adapter 338	TCACAATTTATT
Barcode Adapter 339	GCCAAGACAGAA
Barcode Adapter 340	CTGTGATTTATT
Barcode Adapter 341	ATTGCGGTGAAA
Barcode Adapter 342	GACAGACCCTTT
Barcode Adapter 343	GATATCTAATTT
Barcode Adapter 344	ATCGATAGTAAA
Barcode Adapter 345	TGAGGACTAATT
Barcode Adapter 346	AATCAGGCTTAA
Barcode Adapter 347	CAGAAACTATTT
Barcode Adapter 348	TGCTGTTCTACC
Barcode Adapter 349	GTCCTACTATTT
Barcode Adapter 350	CATCATTCTACC
Barcode Adapter 351	CGACCTGTATTT
Barcode Adapter 352	GCGTGAACCTCC
Barcode Adapter 353	GCTGGTGTATTT
Barcode Adapter 354	CTACACCCTAAA
Barcode Adapter 355	ATCAATGTATTT
Barcode Adapter 356	TATGTAACCTACC
Barcode Adapter 357	TAGTTTGTATTT
Barcode Adapter 358	ATACAAACTCCC
Barcode Adapter 359	GAGCCCATAATT
Barcode Adapter 360	CGCGTTTCTTAA
Barcode Adapter 361	CCTATTACGTCC
Barcode Adapter 362	GCGGAAGTCCTT
Barcode Adapter 363	ACGTTCCGGTCC
Barcode Adapter 364	GGACATGACCGG
Barcode Adapter 365	GGGCTTCAGGCC
Barcode Adapter 366	ACATACGGCCGG
Barcode Adapter 367	CATACACCATCC
Barcode Adapter 368	GGATTGGTTCGG
Barcode Adapter 369	GGCACGTCAGCC

Barcode Adapter 370	ACTGGAATTCGG
Barcode Adapter 371	AAGTCGAGAGCC
Barcode Adapter 372	GGCATCGATCTT
Barcode Adapter 373	GGCAGTCTCTGG
Barcode Adapter 374	CCTGCATCGCAA
Barcode Adapter 375	GGGAGCTAGTGG
Barcode Adapter 376	AACTAGCGCCAA
Barcode Adapter 377	CCCTCCTAGTGG
Barcode Adapter 378	GGTCGGCGCCAA
Barcode Adapter 379	CCGAACCGGTGG
Barcode Adapter 380	GGAGTGTACCAA
Barcode Adapter 381	GGCCAACCACGG
Barcode Adapter 382	ACTTTGGTTGCC
Barcode Adapter 383	AACTGCATGCGG
Barcode Adapter 384	GGGAAGGCCTAA