

NEXTflex™ qRNA-Seq™ v2 Barcoded Primers

NEXTflex™	Forward Sequence	Reverse Complement
Index 1	GTGTTCTA	TAGAACAC
Index 2	AAACATCG	CGATGTTT
Index 3	CATCAAGT	ACTTGATG
Index 4	GTCTGTCA	TGACAGAC
Index 5	AACAACCA	TGGTTGTT
Index 6	CATACCAA	TTGGTATG
Index 7	TGGAACAA	TTGTTCCA
Index 8	ACCTCCAA	TTGGAGGT
Index 9	ACACGACC	GGTCGTGT
Index 10	CTGAGCCA	TGGCTCAG
Index 11	ACACAGAA	TTCTGTGT
Index 12	TATCAGCA	TGCTGATA
Index 13	CTGTAGCC	GGCTACAG
Index 14	TGGCTTCA	TGAAGCCA
Index 15	CCAGTTCA	TGAACTGG
Index 16	AGATGTAC	GTACATCT
Index 17	ACAGCAGA	TCTGCTGT
Index 18	ACAAGCTA	TAGCTTGT
Index 19	ACGCTCGA	TCGAGCGT
Index 20	CACTTCGA	TCGAAGTG
Index 21	AAGACGGA	TCCGTCTT
Index 22	AACGCTTA	TAAGCGTT
Index 23	GAATCTGA	TCAGATTC
Index 24	ATCCTGTA	TACAGGAT
Index 25	AGTACAAG	CTTGTA CT
Index 26	ACCACTGT	ACAGTGGT
Index 27	AACTCACC	GGTGAGTT
Index 28	CGCATACA	TGTATGCG
Index 29	AGGCTAAC	GTTAGCCT
Index 30	GCCAAGAC	GTCTTGGC
Index 31	ATGCCTAA	TTAGGCAT
Index 32	GGTGCGAA	TTCGCACC
Index 33	CAATGGAA	TTCCATTG
Index 34	ATCATTCC	GGAATGAT
Index 35	CAAGACTA	TAGTCTTG
Index 36	AACCGAGA	TCTCGGTT
Index 37	AAGAGATC	GATCTCTT
Index 38	CAAGGAGC	GCTCCTTG
Index 39	CCGTGAGA	TCTCACGG
Index 40	CCGAAGTA	TACTTCGG

Index 41	ACAGATTC	GAATCTGT
Index 42	GACTAGTA	TACTAGTC
Index 43	CGACTGGA	TCCAGTCG
Index 44	CGGATTGC	GCAATCCG
Index 45	ACATTGGC	GCCAATGT
Index 46	AATGTTGC	GCAACATT
Index 47	CAACCACA	TGTGGTTG
Index 48	CCGACAAC	GTTGTCCG
Index 49	AACGTGAT	ATCACGTT
Index 50	GGAGAACA	TGTTCTCC
Index 51	AAGGACAC	GTGTCCTT
Index 52	GTACGCAA	TTGCGTAC
Index 53	AGAGTCAA	TTGACTCT
Index 54	TTCACGCA	TGCGTGAA
Index 55	CACCTTAC	GTAAGGTG
Index 56	GAGCTGAA	TTCAGCTC
Index 57	GCCACATA	TATGTGGC
Index 58	AGCACCTC	GAGGTGCT
Index 59	GCTAACGA	TCGTTAGC
Index 60	AGTCACTA	TAGTGA CT
Index 61	GATGAATC	GATTCATC
Index 62	TGAAGAGA	TCTCTTCA
Index 63	CGCTGATC	GATCAGCG
Index 64	GAGTTAGC	GCTAACTC
Index 65	ATTGGCTC	GAGCCAAT
Index 66	CGAACTTA	TAAGTTCG
Index 67	AATCCGTC	GACGGATT
Index 68	ATAGCGAC	GTCGCTAT
Index 69	TAGGATGA	TCATCCTA
Index 70	CAGCGTTA	TAACGCTG
Index 71	GCTCGGTA	TACCGAGC
Index 72	CGACACAC	GTGTGTCG
Index 73	CAGATCTG	CAGATCTG
Index 74	GATAGACA	TGTCTATC
Index 75	AAGGTACA	TGTACCTT
Index 76	CCTAATCC	GGATTAGG
Index 77	AGATCGCA	TGCGATCT
Index 78	ACGTATCA	TGATACGT
Index 79	AGCAGGAA	TTCCTGCT
Index 80	ACTATGCA	TGCATAGT
Index 81	GCGAGTAA	TFACTCGC
Index 82	AGTGGTCA	TGACCACT
Index 83	CCATCCTC	GAGGATGG
Index 84	CTGGCATA	TATGCCAG
Index 85	AGCCATGC	GCATGGCT
Index 86	CCTCTATC	GATAGAGG
Index 87	TCCGTCTA	TAGACGGA

Index 88	GTCGTAGA	TCTACGAC
Index 89	GAACAGGC	GCCTGTTC
Index 90	CTCAATGA	TCATTGAG
Index 91	CCTCCTGA	TCAGGAGG
Index 92	ATTGAGGA	TCCTCAAT
Index 93	GACAGTGC	GCACTGTC
Index 94	CTAAGGTC	GACCTTAG
Index 95	TGGTGGTA	TACCACCA
Index 96	TCTCACA	TGTGAAGA

©Bioo Scientific Corp. • 2015 v1.1