



BIOO SCIENTIFIC®

THE NGS EXPERTS™

NEXTflex™ qRNA-Seq™ v2 Molecular Indexes

NEXTflex™ Molecular Indexes	Index Sequence
Molecular Index A1	AAC GCC AT
Molecular Index A2	AAG GTA CG
Molecular Index A3	AAT TCC GG
Molecular Index A4	ACA CAG AG
Molecular Index A5	ACA CTC AG
Molecular Index A6	ACA CTG TG
Molecular Index A7	ACA GGA CA
Molecular Index A8	ACC TGT AG
Molecular Index A9	ACG AAG GT
Molecular Index A10	ACG ACT TG
Molecular Index A11	ACG TCA AC
Molecular Index A12	ACG TCA TG
Molecular Index A13	ACT GTC AG
Molecular Index A14	ACT GTG AC
Molecular Index A15	AGA CAC TC
Molecular Index A16	AGA GGA GA
Molecular Index A17	AGC ATC GT
Molecular Index A18	AGC ATG GA
Molecular Index A19	AGC TAC CA
Molecular Index A20	AGC TCT AG
Molecular Index A21	AGG ACA AC
Molecular Index A22	AGG ACA TG
Molecular Index A23	AGG TTG CT
Molecular Index A24	AGT CGA GA
Molecular Index A25	AGT GCT GT
Molecular Index A26	ATA AGC GG
Molecular Index A27	ATC CAT GG
Molecular Index A28	ATC GAA CC
Molecular Index A29	ATC GCG TA
Molecular Index A30	ATC GTT GG
Molecular Index A31	CAA CGA TC
Molecular Index A32	CAA CGT TG
Molecular Index A33	CAA CTG GT
Molecular Index A34	CAA GTC GT
Molecular Index A35	CAC ACA CA
Molecular Index A36	CAG TAC TG
Molecular Index A37	CAT CAG CA
Molecular Index A38	CAT CGT TC
Molecular Index A39	CCA AGG TT
Molecular Index A40	CCT AGC TT

Molecular Index A41	CGA TTA CG
Molecular Index A42	CGC CTA TT
Molecular Index A43	CGT TCC AT
Molecular Index A44	CGT TGG AT
Molecular Index A45	CTA CGT TC
Molecular Index A46	CTA CTC GT
Molecular Index A47	CTA GAG GA
Molecular Index A48	CTA GGA AG
Molecular Index A49	CTA GGT AC
Molecular Index A50	CTC AGT CT
Molecular Index A51	CTG ACT GA
Molecular Index A52	CTG AGT GT
Molecular Index A53	CTG ATG TG
Molecular Index A54	CTG TTC AC
Molecular Index A55	CTT CGT TG
Molecular Index A56	GAA CAG GT
Molecular Index A57	GAA GAC CA
Molecular Index A58	GAA GTG CA
Molecular Index A59	GAC ATG AG
Molecular Index A60	GAG AAG AG
Molecular Index A61	GAG AAG TC
Molecular Index A62	GAT CCT AG
Molecular Index A63	GAT GTC GT
Molecular Index A64	GCC GAT AT
Molecular Index A65	GCC GAT TA
Molecular Index A66	GCG GTA TT
Molecular Index A67	GGA ATT GG
Molecular Index A68	GGA TAA CG
Molecular Index A69	GGC CTA AT
Molecular Index A70	GGC GTA TT
Molecular Index A71	GTC TTG TC
Molecular Index A72	GTG ATG AG
Molecular Index A73	GTG ATG TC
Molecular Index A74	GTG TAC TG
Molecular Index A75	GTG TAG TC
Molecular Index A76	GTT CAC CT
Molecular Index A77	GTT CTG CT
Molecular Index A78	GTT GTC GA
Molecular Index A79	TAC GAA CC
Molecular Index A80	TAG CAA GG
Molecular Index A81	TAG CTA GC
Molecular Index A82	TAG GTT CG
Molecular Index A83	TAT AGC GC
Molecular Index A84	TCA GGA CT
Molecular Index A85	TCC ACA TC
Molecular Index A86	TCG ACT TC
Molecular Index A87	TCG TAG GT

Molecular Index A88	TCG TCA TC
Molecular Index A89	TGA GAC TC
Molecular Index A90	TGA GAG TG
Molecular Index A91	TGA GTG AG
Molecular Index A92	TGC TTG GA
Molecular Index A93	TGG AGT AG
Molecular Index A94	TGT GTG TG
Molecular Index A95	TTC GCC TA
Molecular Index A96	TTC GTT CG